



The Microline® T.F.C.–335 R.O. Drinking Water System will reduce the following:

- Arsenic V
- Barium
- Cadmium
- Chromium III
- Chromium VI
- Copper
- Cyst
- Fluoride
- Lead
- Nitrate
- Nitrite
- Radium 226/228
- Selenium
- TDS

* For specific percentages of reduction, consult the performance data sheet. Do not use with water that is microbiologically unsafe or of unknown quality, without adequate disinfection before or after the system. This system conforms to NSF/ANSI Standard 58 for pentavalent arsenic reduction. See the Performance Data Sheet and Arsenic Facts section for an explanation of reduction performance. Systems certified for cyst reduction may be used on disinfected water that may contain filterable cysts.


This system is acceptable for treatment of influent concentrations of no more than 27 mg/L nitrate and 3 mg/L nitrite in combination measured as N and is certified for nitrate/nitrite reduction only for water supplies with a pressure of 280 kPa (40 psig) or greater.

On a yearly basis it is recommended that you have your water tested by a water treatment professional or an independent laboratory to verify your system is performing satisfactorily.



This system has been Tested and Certified by NSF International against NSF/ANSI Standard 58.

NSF is the nationally recognized agency for the certification of reverse osmosis systems. Products with the NSF mark are certified to pass a series of stringent independent tests. Look for the NSF mark. It is your assurance the system meets its claims.



Your local water treatment professional:

U.S. Patent Numbers: 4,806,912\RE35,252
Canadian Patent Numbers: 2,078,209\2,142,424
Other U.S. and Foreign Patents Pending



4462 Duraform Lane • Windsor • WI 53598-9716
www.microlinero.com

• A Simple, Basic,
• Pure, Necessity... Water.

Because water is vital...


Water is a daily, life-sustaining requirement. It is one of the pure and basic necessities in life.

Ensuring that the water we drink everyday is the safe, natural water that our bodies demand is as vital as the water itself.

Unfortunately, not all of the water we use from our taps is the pure, fresh water that our bodies desire.

The Microline® R.O. Drinking Water System effectively reduces harmful impurities,* turning ordinary tap water into the water that nature intended for us to drink...fresh, natural, life-sustaining water.

Microline® R.O.
Drinking Water Systems



Water.
The crucial
element.



The element of quality.

MICROLINE® R.O. DRINKING WATER SYSTEMS



Model T.F.C.-335

Some of the crystal clear benefits of Microline® R.O. water:

- Delicious, sparkling-clear drinking water.
- Pristine, flavorful coffee, tea and juice.
- Cleanly rinsed fresh fruits and vegetables.
- Crystalline, harder and clearer ice cubes.
- Better tasting soups, sauces and meals when prepared with R.O. water.
- Cost effective. The daily cost of bottled water will no longer be a concern.
- Convenient. Fresh, clean water, ready at your tap.

Microline® R.O. water is water the way that it's supposed to be. Water that looks, smells and tastes like drinking water.

• How Does The • Microline® R.O. System • Work?

The “R.O.” of the Microline® System is the secret.

“R.O.” is **Reverse Osmosis**. It may sound technical, but osmosis is a natural, organic phenomenon, a process that occurs in nature on a continuous basis. Vegetation, like trees, plants and flowers attain their nutrients by using osmosis to draw water from the soil.

The Microline® process of Reverse Osmosis (R.O.) works like this:

The pressure from a household tap forces water through a Microline® semipermeable Membrane. This Membrane separates the everyday household water at the molecular level. The Membrane acts like a filter to isolate impurities,* reducing them so that the water becomes purer, **Reverse Osmosis** water, with substantially reduced dissolved solids and impurities*. This cleaner, more refined water is then stored in a holding tank, ready at your convenience.

• The Microline® Water • Quality Monitor

How do I know if the R.O. System is working properly?

Microline® has taken the guess work out. An optional water quality monitor can be integrated into the system so you can feel more confident that your system is working by simply pushing a button. When you push the water quality monitor button in the cover, a green light shows that your R.O. is working fine. An amber light will show when service may be required.

Systems without the water quality monitor have a water quality test kit that allows you to test your R.O. water with a local water treatment professional or send it out to an independent laboratory for analysis.

*For specific percentages of reduction, consult the performance data sheet.

• The T.F.C.–335 R.O. • Drinking Water System:

The advantage of the T.F.C.–335 R.O. System is that it combines high quality water with flexibility. This system can be installed on well or municipal water systems, to bring R.O. water to the country or urban home. The Thin Film Composite (T.F.C.) reverse osmosis Membrane rejects a higher percentage of impurities* as compared to a C.T.A. Membrane, to bring you the most advanced R.O. technology at an affordable price. Here's how it works:

Water first enters the T.F.C.–335 R.O. System through the Sediment/Carbon Prefilter. This specialized filter protects the automatic shut off and Membrane from clogging with debris, and also reduces chlorine, to protect the refined T.F.C. Membrane.

Water then travels to the operational center of the system—the T.F.C. Membrane. This is where the “Reverse Osmosis” of the Microline® R.O. System happens.

After the Membrane, water is then routed to the holding tank, ready for use. When you turn on the Microline® faucet and draw water, it then goes through its final stage of filtration. The water travels through the Activated Carbon Filter, designed to reduce any remaining tastes and odors before reaching your glass. This carbon imparts the final “sweet” polish to your water.

